


Solution for COVID - 19

eKuore Pro stethoscope

The stethoscope that allows remote auscultation



Use cases



Basic auscultation

It is the most basic form of auscultation, just as you would do with a traditional stethoscope. Using headphones directly connected to the eKuore Pro stethoscope.



Bluetooth auscultation

eKuore Pro enables wireless auscultation using an additional Bluetooth module that sends the auscultation directly to a Bluetooth hearing aids or headphones



Remote auscultation

To reduce the risk of contagion, a sanitary can make the recording and the doctor can listen and see the phonogram simultaneously or later.



Remote monitoring

In cases where the patient is isolated, health professionals can monitor with Personal Protection Equipment (PPE) and the doctor can listen remotely.

What does it do?



REMOTE AUSCULTATION

In the absence of a cable, the medical professional can listen to the auscultations even from an isolated room and without the need to wear a special suit to avoid infections.



REPLACEABLE CHESTPIECES

The chestpiece is the most critical piece of a stethoscope, because it is the part that comes in contact with the infected patient. eKuore Pro has interchangeable and disinfectable chestpieces.



REMOTE MONITORING

The patients can be quickly auscultated from their home and send the sound to healthcare professionals, thus reducing the risk of contagion to the minimum.