



4-year guarantee

# PARI BOY® Classic

## Technical datasheet

**Efficient inhalation system  
for treating diseases of the lower airways  
(chronic and acute)**



**PIF Control  
Patented system**

**Technical data:**

**Alternating current:** 220 – 240V, 50Hz, **Pressure:** 1.6bar, **Weight:** 1.7kg, **Housing dimensions:** 18.5 x 13 x 15 cm (WxHxD), **Fill volume:** min. 2 ml, max. 8 ml, **Noise volume:** 54 db(A)

**The PARI BOY Classic includes:**

- PARI LC SPRINT nebuliser with PIF Control
- Blue nozzle insert for children aged 4+ and for adults
- Mouthpiece with valve system
- Connection tube
- With corresponding accessories, suitable for all age groups



**PARI PIF Control – Patented system:**







PARI nebulisers are equipped with the PIF Control System. The PIF Control System is a breath sensory feedback system for the patient that automatically limits the inhaled air flow (PIF = Peak Inspiratory Flow). If the patient breathes in too quickly the valve reduces the size of the aperture in the nebuliser chamber. The patient perceives this as a resistance – a sign that he should breathe in more slowly for the rest of the session. So, the PIF Control System trains a slow, controlled inspiration resulting in effective drug deposition.

Aerosol characteristics			
Nozzle insert	Total Output Rate in mg/min	Median Mass Diameter in µm	Mass fraction <5µm
Inspiratory flow adults (20l/m)			
Blue	600 mg/min	3.5 µm	67 %

Measurement with Malvern MasterSizer X at 23 °C and 50 % relative humidity. Nebulised medium: 0.9% NaCl (5 ml).

## Nebulisers are available with different nozzle inserts, enabling droplet sizes to be tailored to each age or indication

- Different nozzle inserts produce different aerosol spectra
- Deposition is targetted to the desired region of the airways
- Can be used by a wide age range of patients and for a large variety of indications
- Efficiency and safety of the PARI LC SPRINT® nebuliser has been investigated in 50 scientific studies

Product	Target group/characteristics	Indications	Equipment	REF (Item No.)
 <p>PARI BOY Classic Year Pack</p>	<p>Blue nozzle insert produces droplets for efficient deposition in the central region of the lungs in older children (from 4 years of age) and adults</p> 	<p>Bronchial asthma/ Allergic asthma, Acute and chronic bronchitis, Atelectasis or bronchiectasis, COPD</p>	<p>PARI LC SPRINT® nebuliser (blue nozzle insert) with PIF Control, connection tubing, 2 air filters (one of them for the previous compressor version), mouthpiece</p>	<p>Western Europe 023G3211 IFU languages: de, fr, it, pt, es, nl, en</p> <p>Eastern Europe 023G3216 IFU languages: bg, cs, el, hr, hu, pl, ro, sk, sl, et, lt, lv</p>
 <p>PARI BOY Junior Year Pack</p>	<p>Yellow nozzle insert produces small droplets for small airways in babies and children. Additional red nozzle insert produces finest droplets for smallest airways.</p> 	<p>Bronchial asthma/ Allergic asthma, Acute and chronic bronchitis, Atelectasis or bronchiectasis</p>	<p>PARI LC SPRINT® Junior nebuliser (yellow nozzle insert) with PIF Control, additional red insert, connection tubing, 2 air filters (one of them for the previous compressor version), mouthpiece</p>	<p>Western Europe 023G3311 IFU languages: de, fr, it, pt, es, nl, en</p> <p>Eastern Europe 023G3316 IFU languages: bg, cs, el, hr, hu, pl, ro, sk, sl, et, lt, lv</p>
 <p>PARI BOY Pro Year Pack</p>	<p>Blue nozzle insert produces droplets for efficient deposition in the central region of the lungs. Additional red nozzle insert produces finest droplets for peripheral areas of the lungs.</p> 	<p>Bronchial asthma/ Allergic asthma, Acute and chronic bronchitis, Atelectasis or bronchiectasis, Pulmonary emphysema, COPD</p>	<p>PARI LC SPRINT® nebuliser (blue nozzle insert) with PIF Control, additional red insert, connection tubing, 2 air filters (one of them for the previous compressor version), mouthpiece</p>	<p>Western Europe 023G3011 IFU languages: de, fr, it, pt, es, nl, en</p> <p>Eastern Europe 023G3016 IFU languages: bg, cs, el, hr, hu, pl, ro, sk, sl, et, lt, lv</p>

## With the right accessories for efficient inhalation for all indications and age groups:

Product	REF (Item No.)
BABY mask incl. BABY bend size 0 (for premature infants)	041G0700
BABY mask incl. BABY bend size 1 (for approx. age 0–1 year)	041G0701
BABY mask incl. BABY bend size 2 (for approx. age 1–3 years)	041G0702
BABY mask incl. BABY bend size 3 (for approx. 3 years and older)	041G0703
Child mask soft "Spiggy"	041G0741*
Adult mask soft	078G5000*
Connection tubing, long	041B4587
PARI Filter/Valve Set	041G0500*
PARI Filter Pad, Pack of 30	041B0522
PARI Filter Pad, Pack of 100	041B0523
PARI Filter Pad, Pack of 1000	041B0524
PARI PEP S System	018G4000*
PARI LC Interrupter	022G1000

\* incl. IFU Western Europe de, en, da, el, es, fi, fr, it, nl, no, pt, sv, tr



PARI GmbH  
 Moosstraße 3, 82319 Starnberg, GERMANY  
 E-Mail: [export@pari.de](mailto:export@pari.de)  
 Phone: +49 (0) 8151 279 220  
[www.pari.com](http://www.pari.com)