

UM-211

All-in-one Blood Pressure Monitor



UM-211 is designed for multipurpose use in hospitals. This model is available for dual measurement mode (auscultatory mode with use of the stethoscope and oscillometric mode).

- Dual measurement mode (Auscultatory and Oscillometric)
- Cuff holder and grip for carrying
- Durable and chemical resistant body
- Professional cuff and 5 cuff sizes are available for use
- Large LCD display with backlight
- Rechargeable battery (300 measurements at full charge)
- Adjustable inflation pressure (AUTO/220/250/280)
- Irregular HeartBeat (IHB) indicator
- Room temperature display

UM-211's auscultatory mode is designed for daily and frequent use. Selectable deflation speed (2.5mmHg / sec. or 5.0mmHg / sec.) and quick start for auscultatory mode are features 'Only from A&D'.



Stand-by display (Room temperature)



Measurement display (with backlight)

Multiple Usages



Outpatient



Stored and charged at nurse station



Easy to carry



Measurement at ward

Fits Various Patients

The newly designed professional cuff of the UM Series ranges from SS to LL (5 sizes), 12cm to 50cm for the arm circumference.

| Cuff Size | Arm Size | REF (for UM-102) | REF (for UM-201& UM-211*) |
|-------------|--------------|---------------------|------------------------------|
| Very Small | 12cm to 17cm | CUF-KW-SS | CUF-KS-SS |
| Small Adult | 16cm to 24cm | CUF-KW-SA | CUF-KS-SA |
| Adult | 22cm to 32cm | CUF-KW-A | CUF-KS-A |
| Large Adult | 31cm to 45cm | CUF-KW-LA | CUF-KS-LA |
| Very Large | 41cm to 50cm | CUF-KW-LL | CUF-KS-LL |

